

# sCF 25/8 - starter cell

for evaluation of E/I characteristics with few efforts

- single cell tester 25cm<sup>2</sup>
- graphite flow field, 5 fold serpentine
- electrical heater and bimetallic switch integrated
- gold plated current collectors
- media connection for 6mm tubes
- fixing of test cell with 8 bolts
- operation temperature up to 70°C

optional:

- H<sub>2</sub> CCM 25cm<sup>2</sup> with Pt catalyst 0.5mg/cm<sup>2</sup> on Nafion 212
- GDLs and gaskets

